
Invacare® Kirilos / Atlante

Kirilos, Atlante Electric, Atlante Hydraulic

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This manual MUST be given to the user of the product.
BEFORE using this product, this manual MUST be read and saved
for future reference.

Yes, you can.®

©2020 Invacare Corporation

All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

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Contents

This manual MUST be given to the user of the product.
BEFORE using this product, this manual MUST be read and saved
for future reference.

1 General	4
1.1 Introduction	4
1.1.1 Symbols in this Document	4
1.2 Service Life	4
1.2.1 Additional information	4
1.3 Limitation of Liability	4
1.4 Warranty Information	4
1.5 Compliance	4
1.5.1 Product-specific standards	4
2 Safety	5
2.1 General safety information	5
2.1.1 Pinch points	5
2.2 Safety Information on Accessories	6
2.3 Safety information on electromagnetic interference	6
2.4 Labels and symbols on the product	6
2.4.1 Identification label	6
2.4.2 Other Labels and Symbols	7
3 Product Overview	8
3.1 Intended use	8
3.2 Main parts of the lift	8
3.3 Accessories	8
4 Setup	9
4.1 General Safety Information	9
4.2 Scope of delivery	9
4.3 Assembling the Lift	9
4.4 Disassembling the Lift	10
5 Usage	11
5.1 General safety information	11
5.2 Locking/Unlocking the rear castors	11
5.3 Raising/Lowering an Electric Lift	11
5.4 Raising/Lowering a Hydraulic Lift	11
5.5 Closing/Opening Legs	11
5.5.1 Closing/Opening legs manually	12
5.6 Emergency functions	12
5.6.1 Performing an Emergency Stop	12
5.6.2 Activating an Emergency Lowering	12
5.7 Charging the battery	12
6 Patient Transfer	14
6.1 General safety information	14
6.2 Preparing to lift	14
6.2.1 Attaching the sling to the lift	15
6.3 Transferring a patient from a bed	15
6.4 Transferring a patient to a bed	16
6.5 Transferring a patient from a wheelchair	16
6.6 Transferring a patient to a wheelchair	17
6.7 Transferring a patient to and from a commode	17
6.8 Lifting a patient from the floor	18
7 Maintenance	19
7.1 General maintenance information	19
7.2 Daily Inspections	19
7.3 Cleaning and Disinfection	19
7.3.1 General safety information	19
7.3.2 Cleaning Intervals	19
7.3.3 Cleaning instructions	19
7.3.4 Disinfection Instructions	19
7.4 Service interval	20
7.4.1 Safety inspection checklist	20
7.5 Lubricating the Lift	20
7.6 Checking and tightening bolts	20
7.7 Maintaining the Hydraulic Pump	21
8 After Use	22
8.1 Transport and Storage	22
8.2 Disposal	22
8.3 Reconditioning	22
9 Troubleshooting	23
9.1 Identifying and repairing faults	23
10 Technical Data	24
10.1 Dimensions and weight	24
10.2 Electrical system	24
10.3 Hydraulic system	25
10.4 Environmental conditions	25
10.5 Materials	25
11 Electromagnetic compatibility (EMC)	26
11.1 General EMC information	26
11.2 Electromagnetic emission	26
11.3 Electromagnetic Immunity	26
11.4 Test specifications for immunity to RF wireless communications equipment	27

1 General

1.1 Introduction

This user manual contains important information about the handling of the product. To ensure safety when using the product, read the user manual carefully and follow the safety instructions.

Note that there may be sections in this document, which are not relevant to your product, since this document applies to all available models (on the date of printing). If not otherwise stated, each section in this document refers to all models of the product.

The models and configurations available in your country can be found in the country-specific sales documents.

Invacare reserves the right to alter product specifications without further notice.

Before reading this document, make sure you have the latest version. You find the latest version as a PDF on the Invacare website.

If you find that the font size in the printed document is difficult to read, you can download the PDF version from the website. The PDF can then be scaled on screen to a font size that is more comfortable for you.

For more information about the product, for example product safety notices and product recalls, contact your Invacare distributor. See addresses at the end of this document.

In case of a serious incident with the product, you should inform the manufacturer and the competent authority in your country.

1.1.1 Symbols in this Document

Symbols and signal words are used in this document and apply to hazards or unsafe practices which could result in personal injury or property damage. See the information below for definitions of the signal words.

WARNING

Indicates a hazardous situation that could result in serious injury or death if it is not avoided.

CAUTION

Indicates a hazardous situation that could result in minor or slight injury if it is not avoided.

IMPORTANT

Indicates a hazardous situation that could result in damage to property if it is not avoided.

Tips and Recommendations

Gives useful tips, recommendations and information for efficient, trouble-free use.

1.2 Service Life

The expected service life of this product is eight years when used daily and in accordance with the safety instructions,

maintenance intervals and correct use, stated in this manual. The effective service life can vary according to frequency and intensity of use.

1.2.1 Additional information

The expected service life is based on an estimated average of 4 lifting cycles per day.

1.3 Limitation of Liability

Invacare accepts no liability for damage arising from:

- Non-compliance with the user manual
- Incorrect use
- Natural wear and tear
- Incorrect assembly or set-up by the purchaser or a third party
- Technical modifications
- Unauthorised modifications and/or use of unsuitable spare parts

1.4 Warranty Information

We provide a manufacturer's warranty for the product in accordance with our General Terms and Conditions of Business in the respective countries.

Warranty claims can only be made through the provider from whom the product was obtained.

1.5 Compliance

Quality is fundamental to the company's operation, working within the disciplines of ISO 13485.

This product features the CE mark, in compliance with the Medical Device Regulation 2017/745 Class I. The launch date of this product is stated in the CE declaration of conformity.

We are continuously working towards ensuring that the company's impact on the environment, locally and globally, is reduced to a minimum.

We only use REACH compliant materials and components.

We comply with the current environmental legislations WEEE and RoHS.

1.5.1 Product-specific standards

The product has been tested and conforms to ISO 10535 (Hoists for the transfer of disabled persons) and all related standards.

For further information about local standards and regulations, contact your local Invacare representative. See addresses at the end of this document.

2 Safety

2.1 General safety information

This section of the manual contains general safety information about your product. For specific safety information, refer to the appropriate section of the manual and procedures within that section.

WARNING!

Risk of injury or damage

- Do not use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material such as user manuals or instruction sheets supplied with this product or optional equipment. If you are unable to understand the warnings, cautions or instructions, contact a healthcare professional, Invacare provider or qualified technician before attempting to use this product.
- Do not make any unauthorized alterations or modifications to the product.

WARNING!

The maximum safe working load must not be exceeded

- Do not exceed the maximum safe working load of this product or used accessories like slings, spreader bars etc. See documentation or the labeling for the stated maximum safe working load.
- The component with the lowest load limit determines the maximum safe working load of the entire system.

WARNING!

Risk of injury or damage

- Improper use of this product may cause injury or damage.
- Do not attempt any transfer without approval of the patient's healthcare professional.
 - Read the instructions in this user manual and observe trained personnel performing transfer procedures. Then practice transfers under supervision and with a capable person acting as a patient.
 - Special care must be taken with people with disabilities who cannot cooperate while being transferred.
 - Do not use the lift as a transport device. It is intended to transfer an individual from one resting surface to another.

WARNING!

Risk of injury or damage

- Inappropriate handling of cables can cause electrical shock and product failure.
- Do not kink, shear or otherwise damage the cables of the product.
 - Make sure that no cables are jammed or damaged, when the product is used.
 - Make sure the wiring is correct and connections are proper.
 - Do not use unauthorized equipment.

WARNING!

Risk of injury or damage

- Excessive moisture will damage the product and may cause electrical shock.
- The patient lift can be used in a bath or shower area but must NOT be used under the shower. The patient must be transferred to a shower chair or use other means for showering.
 - If the patient lift is used in a moisture environment, ensure that the patient lift is wiped clean of any moisture after use.
 - Do not plug or unplug the power cable in a moisture environment or with wet hands.
 - Do not store the product in a damp area or damp condition.
 - Periodically inspect all components of the product for signs of corrosion or damage. Replace parts that are corroded or damaged.
 - See 10.4 *Environmental conditions*, page 25.

WARNING!

Risk of injury or damage

- Ignition sources can cause burns or fire.
- Patient transfer must be performed with a safety space between the lift and possible ignition sources (heater, stove, fireplace, etc.)
 - The patient and assistants must not smoke during transfer.
 - The sling must not be placed over heat sources (heater, stove, fireplace, etc.)

WARNING!

Risk of injury or damage

- To avoid injury or damage when operating the product:
- Close supervision is necessary when the product is used near children or pets.
 - Do not let children play with the product.

CAUTION!

Risk of injury or damage


- The product can get hot when exposed to sunlight or other heat sources.
- Do not expose the product to direct sunlight for prolonged periods.
 - Keep the product away from heat sources.

IMPORTANT!

- Accumulation of lint, dust and other dirt can impair the product.
- Keep the product clean.

2.1.1 Pinch points

WARNING!

Risk of Injury

Pinch points are present in several locations on the lift and fingers could be pinched.
 – Always keep hands and fingers clear of moving parts.

2.2 Safety Information on Accessories

CAUTION!

Risk of Injury

Non-original or incorrect accessories may affect the function and safety of this product.
 – Due to regional differences, refer to your local Invacare website or catalog for available accessories or contact your Invacare provider.
 – See the manual delivered with the accessory for further information and instructions.
 – Use only original accessories for the product in use. Under certain circumstances, the use of slings from other manufacturers is possible. See additional information in this section.

CAUTION!

Compatibility of Slings with Attachment System

Invacare uses a common attachment system based on hooks and loops. Loops on the slings are attached to hooks on spreader bars (coat hanger style). Therefore, suitable slings of other manufacturers can be used on this lift as well.
 – Use only slings with loop attachments that are suitable for spreader bars (coat hanger style) with hooks.
 – Do not use slings designed for "Keyhole or Clip Attachment Systems" or "Tilting Frame Systems".

To select the appropriate sling, a risk assessment must be performed by a healthcare professional. The risk assessment must consider:

- The patient's weight, size, physical ability and medical condition.
- The type of transfer and the environment.
- The compatibility to the other lifting equipment used.

2.3 Safety information on electromagnetic interference

WARNING!

Risk of malfunction due to electromagnetic interference

Electromagnetic interference between this product and other electrical equipment can occur and disturb the electrical adjustment functions of this product. To prevent, reduce or eliminate such electromagnetic interference:

- Only use original cables, accessories and spare parts, to not increase electromagnetic emission or reduce electromagnetic immunity of this product.
- Do not use portable radio frequency (RF) communications equipment closer than 30 cm to any part of this product (including cables).
- Do not use this product near active high-frequency surgical equipment and the RF shielded room of a system for magnetic resonance imaging, where the intensity of electromagnetic disturbances is high.
- If disturbances occur, increase the distance between this product and the other equipment or switch it off.
- Refer to the detailed information and follow the guidance in chapter 11 *Electromagnetic compatibility (EMC)*, page 26.

WARNING!

Risk of malfunction

Electromagnetic interference may cause improper operation.


- Do not use this product adjacent to or stacked with other electrical equipment. If such use is necessary this product and the other equipment must be closely observed to verify that they are operating normally.


2.4 Labels and symbols on the product

2.4.1 Identification label

The identification label contains the main product information, including technical data.

SN	Serial number
REF	Reference number
	Manufacturer

	Date of manufacture
	Max. safe working load
	Class II equipment (Kirilos and Atlante Electric only)
	Type B applied part (Kirilos and Atlante Electric only)
	WEEE conform (Kirilos and Atlante Electric only)
	Read user manual (Atlante Hydraulic only)
	European Conformity
	Medical device

	
	Total weight of the product with the maximum safe working load applied
	Temperature limit
	Humidity limitation
	Atmospheric pressure limitation
	Transport and storage conditions
	Operation conditions

Abbreviations for technical data:

- I_{in} = Incoming Current
- U_{in} = Incoming Voltage
- Int. = Intermittence
- AC = Alternating Current
- Max = maximum
- min = minute

For more information about technical data, see *10 Technical Data, page 24*.

2.4.2 Other Labels and Symbols

(Kirilos and Atlante Electric only)

Read the user manual before using this product and follow all instructions for safety and use.

Do not apply any side forces to the actuator (e.g. by using it as a push bar) as this might result in damage and malfunction.

Do not remove this label and manipulate the underlying screw as this might impair the safety of the product
The warranty of the product is voided if this label is removed or broken

3 Product Overview

3.1 Intended use

The mobile patient lift is a transfer device and is intended to transfer and position an individual from one resting surface to another. For example:

- Between the bed and a wheelchair
- To and from the toilet
- Lowering and raising patients to/from the floor

The maximum safe working load is stated in *10 Technical Data, page 24*

The mobile patient lift is designed to be used indoors on a level surface, in hospitals, nursing facilities and domestic areas. The mobile patient lift can be turned (rotated) in place for transfers with limited floor space.

A healthcare professional or private person who has received proper training is the intended operator of this product.


Indications

The mobile patient lift is indicated for completely or partially immobile patients, who cannot be transferred with other types of lifts or transfer aids.

All position changes are possible without assistance of the patient.

There are no known contraindications for this product.

3.2 Main parts of the lift

Ⓐ	Base
Ⓑ	Leg
Ⓒ	Castor
Ⓓ	Rear castor with brake
Ⓔ	Boom
Ⓕ	Mast
Ⓖ	Spreader bar
Ⓗ	Hook for sling
Ⓘ	Actuator (Kirilos and Atlante Electric only)
Ⓝ	Control unit with battery (Kirilos and Atlante Electric only)
Ⓚ	Emergency Stop (Kirilos and Atlante Electric only)
Ⓛ	Push bar
Ⓜ	Hand control (Kirilos and Atlante Electric only)
Ⓝ	Leg spreader handle
Ⓞ	Hydraulic pump – replaces actuator Ⓘ for Atlante Hydraulic
Ⓟ	Manual emergency lowering device (Kirilos and Atlante Electric only)

3.3 Accessories

- ⓘ Due to regional differences, refer to your local Invacare website or catalog for details about available accessories or contact your Invacare provider.

Sling models with loop attachments, suitable for spreader bars (coat hanger style) with hooks:

- Full body support slings - without head support
- Full body support slings - with head support
- Slings for dress/ toileting - with or without head support
- Slings for amputee

4 Setup

4.1 General Safety Information

WARNING!
Risk of Injury or Damage

Damaged parts may affect the safety of the product.

- Check all parts for shipping damage before use.
- In case of damage, do not use the equipment and contact your Invacare provider for further instructions.

WARNING!
Risk of Injury or Damage


Improper assembly may cause injury or damage.

- Assembly must be performed only by a qualified technician.
- Use only Invacare parts in the assembly of this patient lift.
- After each assembly, check that all fittings are properly tightened and all parts have the correct function
- Do not overtighten the mounting hardware. This may damage the mounting brackets.

If there are any issues or questions during assembly, contact your Invacare provider.

4.2 Scope of delivery

The items listed in the table are included in your package.

(A)	Mast assembly including boom, spreader bar and actuator or hydraulic pump (1 piece)
(B)	Battery (1 piece, Kirilos and Atlante Electric only)
(C)	Hand control (1 piece, Kirilos and Atlante Electric only)
(D)	Leg spreader handle (1 piece)
(E)	Base with legs (1 piece)

(F) - not shown	Cable for charging (1 piece, Kirilos and Atlante Electric only)
(G) - not shown	Sling (1 piece)
(H) - not shown	User Manual (1 piece)


4.3 Assembling the Lift

WARNING!

- Check all parts for visible defects or damage before assembly. In case of any damage, do not use the product and contact Invacare.
- Make sure the emergency stop is activated before assembly or disassembly.
- Take care when lifting components during assembly. Some parts are heavy. Always remember to adopt the correct lifting position.

Perform unpacking and assembly operation at floor level.

Assembling the mast to the base


1. Put the base (A) on a level surface and lock the rear castors (B).
2. Loosen the mast knob (C), but leave it screwed into the base (A).
3. Insert the tube of the mast (D) into the hole of the base (A).
4. While supporting the mast assembly (D), tighten the mast knob (C).
5. Place the leg spreader handle (E) onto its support (F) and push the spring pin, to lock the handle (E) on its support (F).

Installing the actuator or hydraulic pump to the boom


2 x 17 mm wrench

1. Loosen the nut **A** and remove bolt **B** from the boom mounting bracket **C**.
2. Place the actuator **D** or hydraulic pump **E** in the boom mounting bracket **C** and align the holes.
3. Reinstall the bolt **B** and tighten the nut **A**.

Wiring

(Kirilos and Atlante Electric only)

1. Plug the hand control cable on socket **A** and motor cable on socket **B** on the control unit.
2. Place the battery **C** on the control unit. Push the top of the battery against the mounting bracket until there is an audible click.

4.4 Disassembling the Lift

1. Put the base on a level surface so that all castors make contact with the floor.
2. Lock the rear castors.
3. Reverse the relevant steps of the assembling procedure for your model.

5 Usage

5.1 General safety information

WARNING!

Risk of injury or damage

Before using the lift with a patient, refer to the following safety information and instructions:

- 2 Safety, page 5
- 6 Patient Transfer, page 14

5.2 Locking/Unlocking the rear castors


WARNING!

Risk of Injury or Damage

The lift could tip and endanger the patient and assistants.


- Invacare recommends that the rear castors be left unlocked during lifting procedures to allow the patient lift to stabilize itself when the patient is initially lifted from a chair, bed or any stationary object.
- Invacare recommends locking the rear castors of the lift only when positioning or removing the sling from around the patient.

- To lock the castor, push down pedal **A** with the foot.
- To unlock the castor, push up pedal **A** with the foot.

5.3 Raising/Lowering an Electric Lift

The hand control is used to raise or lower the lift.


1. To raise the lift — Press and hold the UP **A** button to raise the boom and the patient.
2. To lower the lift — Press and hold the DOWN **B** button to lower the boom and the patient.

Release the button to stop raising or lowering the lift.

5.4 Raising/Lowering a Hydraulic Lift

The pump handle and control valve are used to raise or lower the hydraulic lift.

Refer to the safety information in the Introduction of the Usage section before performing this procedure.

There are two controls on the pump assembly:

- The control valve **A**
- The pump handle **B**

To raise the lift:

1. Turn the control valve **A** clockwise to the CLOSED position.
2. Move the pump handle **B** up and down to elevate the boom **C**.

To lower the lift:

1. Turn the control valve **A** counterclockwise to the OPEN position.
2. The weight of the patient in the sling **D** will cause the boom to descend.
3. If there is no patient in the sling, pull down on the boom **C**.

The rate of descent of the boom is controlled by the amount that the control valve is open.

A safety gate is part of the hydraulic system that controls the maximum descent of the boom regardless of how far the control valve is opened.

5.5 Closing/Opening Legs


WARNING!

Risk of Injury

The lift could tip and endanger the patient and assistants.

- The legs of the lift must be in the maximum open position for optimum stability and safety. If it is necessary to close the legs of the lift to maneuver the lift under a bed, close the legs of the lift only as long as it takes to position the lift over the patient and lift the patient off the surface of the bed. When the legs of the lift are no longer under the bed, return the legs of the lift to the maximum open position.

5.5.1 Closing/Opening legs manually

1. To open the legs, pull the leg spreader handle **A** off and then to the left until it locks in the notch of the bracket.
2. To close the legs, push the leg spreader handle **A** off and then to the right until it locks in the notch of the bracket.

 The shifter handle **MUST** lock into its mounting slot to lock the legs in the full closed position.

WARNING!
Risk of Injury


The lift could tip and endanger the patient and assistants.

- If the leg spreader handle is **NOT** positioned completely into its mounting slot, **DO NOT** use the patient lift until the leg spreader handle is properly seated and the legs of the patient lift are locked in place.

5.6 Emergency functions

5.6.1 Performing an Emergency Stop

(Kirilos and Atlante Electric only)


1. Press the red emergency button **A** on the control unit to stop the boom and patient from raising or lowering.
2. To reset, rotate the emergency button clockwise.

5.6.2 Activating an Emergency Lowering

Activating an Emergency Lowering on the Control Unit

(Kirilos and Atlante Electric only)

If the hand control fails, the boom can be lowered by using the circular switch for emergency release. Use a sharp object, such as a pencil to press the button.

1. Lower the boom by pressing and holding button **A** at the front of the control unit.
2. Stop lowering the boom by releasing the button.

Activating a Mechanical Emergency Lowering

(Kirilos and Atlante Electric only)


If the emergency lowering on the control unit is not working, a mechanical emergency lowering is available as a back-up. This may happen in case of partial or total power failure, or if the battery runs down during use.

WARNING!
Risk of Injury or Damage

Repeated mechanical emergency lowering leads to a high increase of the lowering speed.

- Only use the mechanical emergency lowering as a back-up when the emergency lowering on the control unit is not working.
- Do not use the mechanical emergency lowering repeatedly in close succession.
- Always lower the patient onto a supporting surface like a bed or a chair.
- After a mechanical emergency lowering was required, the lift must be checked to resolve the malfunction of the control unit.


1. Locate the emergency release knob **A** at the top of the actuator piston.
2. Rotate the emergency release knob **A** clockwise.

The boom will lower slowly without additional downward force on the boom.

Activating an Emergency Lowering on the Hydraulic Pump

(Atlante Hydraulic only)

1. Turn the control valve **A** counterclockwise to the OPEN position.

5.7 Charging the battery

(Kirilos and Atlante Electric only)


! **IMPORTANT!**

- Make sure the emergency stop is not activated while charging the battery.
- Make sure that charging takes place in a room with good air ventilation.
- Electrical functions are not working when the lift is connected to power supply.
- Do not attempt to use the lift if the battery housing is damaged.
- Replace a damaged battery housing before further use.
- Do not move the lift without unplugging from the power outlet.

It is recommended to charge the battery daily to ensure optimal use of the lift and prolong the life of the battery. Furthermore, it is recommended to charge the battery before first use.

Control unit

The control unit is equipped with a sound signal. A beep indicates that the battery has low capacity, but lowering the patient is still possible. It is recommended to charge the batteries as soon as the sound signal is heard.

1. Plug the cable for charging © into a power outlet.
 - i** The battery will charge in approximately 4 hours. The charger stops automatically when the batteries are fully charged. The right yellow charge diode (A) will continuously light during charging and switch off when fully charged. The left green diode (B) will light continuously while the control unit is connected to the mains, and lights up when any button on the hand control is pressed.
2. Disconnect the power cord from the power outlet after the battery has been fully charged.

6 Patient Transfer

6.1 General safety information

WARNING!

Risk of Injury or Damage

Improper use of this product may cause injury or damage.

- Before transferring to a stationary object (wheelchair, bed, commode or other surface) check that the weight capacity can withstand the patient's weight.
- If applicable the wheel or castor locks of the stationary object (wheelchair, bed etc.) must be in a locked position before lowering the patient onto or lifting the patient off.

WARNING!

Risk of Injury or Damage

Improper use of this product may cause injury or damage.

- Use the push bars on the mast at all times to push or pull the lift.
- Avoid using the lift on an incline. Invacare recommends that the product only be used on a flat surface.
- During transfer, with the patient suspended in a sling attached to the lift, DO NOT roll lift over uneven surfaces that could cause the lift to tip over.

WARNING!

Risk of Injury or Damage

Damage to parts of the lift (hand control, castors, etc.) caused by impact with the floor, walls or other stationary objects may cause damage to the product and lead to injury.

- DO NOT allow parts of the lift to impact the floor, walls or other stationary objects.
- ALWAYS store the hand control properly when not in use.

WARNING!

Risk of Injury

The spreader bar can move suddenly and cause injury.

- When positioning lift, be aware of the position of the spreader bar and the patient.

WARNING!

Risk of Entrapment or Strangulation

The hand control cord can cause injury if improperly positioned and secured.

- ALWAYS be aware of the location of the hand control cord in relation to the patient and caregivers.
- DO NOT allow the hand control cord to become entangled around the patient and caregivers.
- The hand control must be secured properly. ALWAYS store the hand control properly when not in use.

WARNING!

Risk of Entrapment or Strangulation

Items in the patient's surroundings can cause entrapment strangulation during lifting. To avoid entrapment or strangulation:

- Before lifting, check that the patient is completely free of his/her surroundings.

WARNING!

Risk of Entrapment

There is a risk of entrapment between the spreader bar hooks and the sling.

- Use caution when lifting.
- NEVER put hands or fingers on or near the hooks when lifting.
- Ensure the patient's hands and fingers are away from the hooks before lifting.

IMPORTANT!

All transfer procedures described in the following can be performed by one (1) assistant. However, Invacare recommends that the procedures are performed by two (2) assistants whenever possible.

6.2 Preparing to lift

WARNING!

Risk of Injury

During transfers and lift operation, the boom can impact the patient or caregivers and cause injury.

- ALWAYS be aware of the position of the boom during transfers.
- Ensure the boom is positioned in a way that it cannot impact the patient or bystanders.
- ALWAYS be aware of your body position in relation to the boom during transfers.

1. Before proceeding, review the information and observe all warnings indicated in 2 *Safety, page 5* and 6 *Patient Transfer, page 14*.
2. Position the patient onto the sling. Refer to your sling user manual.
3. Unlock the rear castors. See 5 *Usage, page 11*.
4. Open the legs of the lift. See 5 *Usage, page 11*.
5. Use the push bars to move the patient lift into position.

WARNING!

– When using the lift in conjunction with beds or wheelchairs, be aware of the position of the lift in relationship to those other devices so that the lift does not become entangled.

- Before positioning the legs of the patient lift under a bed, make sure that the area is clear of any obstructions.

6. Lower the patient lift for easy attachment of the sling.
7. Lock the rear castors. See 5 *Usage, page 11*
8. Attach the sling. See 6.2.1 *Attaching the sling to the lift, page 15*.

6.2.1 Attaching the sling to the lift

WARNING!
Risk of Injury


- Usage of wrong or damaged slings can cause the patient to fall or cause injury to assistants.
- Use an Invacare approved sling that is recommended by the individual's doctor, nurse or medical assistant for the comfort and safety of the individual being lifted.
 - Invacare slings and patient lift accessories are specifically designed to be used in conjunction with Invacare patient lifts.
 - After each laundering (in accordance with instructions on the sling), inspect sling(s) for wear, tear, and loose stitching.
 - Bleached, torn, cut, frayed, or damaged slings are unsafe and could result in injury. Discard immediately.
 - Do not alter slings.

WARNING!
Risk of Injury

- Improperly attached or adjusted slings can cause the patient to fall or cause injury to assistants.
- Be sure to check the sling attachments each time the sling is removed and replaced, to ensure that it is properly attached before the patient is removed from a stationary object (bed, chair or commode).
 - Do not use any kind of plastic back incontinence pad or seating cushion between patient and sling material that may cause the patient to slide out of the sling during transfer.
 - Make sure that there is sufficient head support when lifting a patient.
 - Position the patient in the sling as directed by the instructions provided with the sling.
 - Adjustments for patient safety and comfort should be made before moving the patient.


The straps of the sling may be equipped with color coded loops providing different lengths to place the patient into different positions. Shorter straps at the shoulders will produce a more vertical lift, which will assist with positioning into a chair or wheelchair. By lengthening the straps at the shoulders, a more reclined position can be achieved which is more suitable for a chair to bed transfer. Always match the loop colors of the corresponding straps on each side of the sling for an even lift of the patient.

1. Place the desired loop of strap (A) over hook (C) on the spreader bar (D).
2. Pull down the strap (A) until the loop is fully seated at the bottom of hook (C).
3. Repeat the steps for each of the remaining straps (A) of the sling (B).


6.3 Transferring a patient from a bed

1. Prepare to lift. See 6.2 *Preparing to lift*, page 14.
2. Attach the sling to the lift. See *Attaching the sling to the lift*.
3. Unlock the rear castors.
- 4.

Lower the bed to the lowest position.


5.

Lift the patient high enough to clear the stationary object with their weight fully supported by the lift.


The boom will stay in position until the DOWN button is pressed.

6.

Before moving the patient, check again to make sure that the sling is properly attached to the hooks of the spreader bar. See *Attaching the sling to the lift*. If any attachments are not properly in place, lower the patient back onto the stationary object and correct the problem.

7.


Using the push bars, move the lift away from the stationary object.

8.

Using the handles on the sling, turn the patient so that he/she faces the assistant operating the patient lift.


9.

Lower the patient so that his/her feet rest on the base of the lift, straddling the mast.

 The lower center of gravity provides stability making the patient feel more secure and the lift easier to move.

10.

Move the patient lift with both hands firmly on the push bars.

11. Proceed with the transfer to another resting surface or stationary object.

6.4 Transferring a patient to a bed

1. Perform the following steps in addition to those required to lift the patient of another resting surface or stationary object
- 2.


Raise or lower the lift to position the patient over the stationary surface. Be sure to raise or lower the patient enough to clear the sides of the stationary object.

3.


Lower the patient onto the stationary surface.

4. Ensure the patient is fully supported by the surface to which you are transferring.
5. Lock the rear castors.
6. Raise the bed to a good working height (usually hip height of caregivers).
- 7.

Detach the sling from the spreader bar.


8. Unlock the rear castors.
- 9.

Move the lift away from the area.

6.5 Transferring a patient from a wheelchair

1. Engage the parking brakes of the wheelchair to prevent movement of the wheelchair.
2. Prepare to lift. See 6.2 *Preparing to lift*, page 14.
3. Ensure the wheelchair is placed between the opened legs of the lift with the patient facing the mast.
4. Attach the sling to the lift. See *Attaching the sling to the lift*.
5. Unlock the rear castors.
- 6.

Raise the boom until there is tension on the straps of the sling and make sure that the sling is properly attached to the hooks of the spreader bar. See *Attaching the sling to the lift*. If any attachments are not properly in place, lower the patient back onto the wheelchair and correct the problem.

7.

Lift the patient high enough to clear the wheelchair surface with their weight fully supported by the lift.

8. Using the push bars, move the lift away from the wheelchair.
9. Proceed with the transfer to another resting surface or stationary object.


6.6 Transferring a patient to a wheelchair

WARNING!
Risk of Injury

- Before transferring, check that the wheelchair weight capacity can withstand the patient’s weight.
- The wheelchair wheel locks **MUST** be in a locked position before lowering the patient into the wheelchair for transport.

1. Perform the following steps in addition to those required to lift the patient of another resting surface or stationary object
2. Engage the parking brakes of the wheelchair to prevent movement of the wheelchair.
3. Position the patient over the seat with their back against the back of the chair.
- 4.

Begin to lower the patient and use the grab handle (on selected models) or sides of the sling to guide the patient’s hips as far back as possible into the seat for proper positioning. This will maintain a good center of balance and prevent the chair from tipping forward.

6.7 Transferring a patient to and from a commode

The Invacare patient lift is **NOT** intended as a transport device. If the bathroom facilities are **NOT** near the bed or if the patient lift cannot be easily maneuvered towards the commode, then the patient **MUST** be transferred to a wheelchair and transported to the bathroom facilities before using the patient lift again to position the patient on a standard commode.


The slings with commode openings are designed to be used with either a commode chair or standard commode.


Before transferring the patient, the patient lift should be guided to the bathroom facilities to check that it can be easily maneuvered towards the commode.

1. Perform the following steps in addition to those required to lift the patient of another resting surface or stationary object
- 2.


Elevate the patient high enough to clear the commode chair arms and have their weight supported by the patient lift. See 5 Usage, page 11

3.

Position the lift so the legs are outside of the commode legs and the push bars of the lift are opposite of the commode and guide the patient over the commode.

4.

Lower the patient onto the commode, leaving the sling attached to the spreader bar hooks. Invacare recommends that the sling remains connected to the spreader bar hooks during the patient’s use of either the commode chair or standard commode.

5.

When complete, recheck for correct sling attachment.

6. Raise the patient off of the commode.
- 7.

When the patient is clear of the commode surface, use the push bars to move the lift away from the commode.

8. Proceed with the transfer to another resting surface or stationary object.


6.8 Lifting a patient from the floor

WARNING! Risk of Injury


The lift can injure a patient who is on the floor or an assistant working on the floor.
– Always take care when moving parts of the lift close to a person on the floor.

1. Determine if the patient has suffered any injuries from a fall. If no medical attention is needed, proceed with the transfer.
- 2.


Position the sling **A** under the patient. Refer to the sling user manual for more information about positioning slings.

3.


Have the patient bend his/her knees and raise his/her head off of the floor supporting the patient's head with a pillow **B**.

4. Unlock the rear castors. See 5 Usage, page 11
5. Open the legs of the lift. See 5 Usage, page 11
- 6.


Position the lift with one leg under the patient's bent knees and the other leg under the patient's head. Keep the sling straps inside of the legs of the lift.

7.

Lower the boom so the spreader bar is directly over the patient's chest.

8.

Attach the sling. See *Attaching the sling to the lift*.

9. Raise the boom to lift the patient off the floor.
10. Proceed with the transfer to another resting surface or stationary object.

7 Maintenance

7.1 General maintenance information

CAUTION!

Risk of injury or damage

- Do not perform any maintenance or service procedures while the product is in use.

Follow the maintenance procedures described in this manual to keep your product in continuous operation.

7.2 Daily Inspections

WARNING!

Risk of injury or damage

Damaged or worn parts can affect the safety of the lift.

- The patient lift should be checked each time it is used.
- Do not use the lift if damage is found or you question the safety of any part of the lift. Contact your Invacare provider immediately and ensure the lift is not used until repairs are made.

Daily Inspection Checklist

- Visually inspect the patient lift. Check all parts for external damage or wear.
- Check all hardware, attachment points and stressed parts, such as slings, hooks and any pivot points for signs of wear, cracking, fraying, deformation or deterioration.
- Verify that the hand control or hydraulic pump is functional.

Kirilos and Atlante Electric only:

- Charge the battery every day the lift is used.
- Check the emergency stop function.

7.3 Cleaning and Disinfection

7.3.1 General safety information

CAUTION!

Risk of contamination

- Take precautions for yourself and use appropriate protective equipment.

CAUTION!

Risk of electric shock and product damage

- Switch off the device and disconnect from mains, if applicable.
- When cleaning electronic components consider their protection class regarding water ingress.
- Make sure that no water splashes to the plug or the wall outlet.
- Do not touch the power socket with wet hands.

IMPORTANT!

Wrong fluids or methods can harm or damage the product.

- All cleaning agents and disinfectants used must be effective, compatible with one another and must protect the materials they are used to clean.
- Never use corrosive fluids (alkalines, acid etc.) or abrasive cleaning agents. We recommend an ordinary household cleaning agent such as dishwashing liquid, if not specified otherwise in the cleaning instructions.
- Never use a solvent (cellulose thinner, acetone etc.) that changes the structure of the plastic or dissolves the attached labels.
- Always make sure that the product is completely dried before taking into use again.

For cleaning and disinfection in clinical or long-term care environments, follow your in-house procedures.

7.3.2 Cleaning Intervals

IMPORTANT!

Regular cleaning and disinfection enhances smooth operation, increases the service life and prevents contamination.

Clean and disinfect the product

- regularly while it is in use,
- before and after any service procedure,
- when it has been in contact with any body fluids,
- before using it for a new user.

7.3.3 Cleaning instructions

IMPORTANT!

- The product does not tolerate cleaning in automatic washing systems, with high-pressure cleaning equipment or steam.

Cleaning the lift

Method: Wipe off with a wet cloth or soft brush.

Max. temperature: 40 °C

Solvent/chemicals: Ordinary household cleaning agent and water.

Drying: Wipe dry with a soft cloth.

Cleaning the sling

Refer to the washing instructions on the sling and to the sling manual for cleaning details.

7.3.4 Disinfection Instructions

Information on recommended disinfectants and methods can be found on <https://vah-online.de/en/for-users>.

In Domestic Care

Method: Follow the application notes for the used disinfectant and wipe-disinfect all accessible surfaces.

Disinfectant: Ordinary household disinfectant.

Drying: Allow the product to air-dry.

In Clinical and Long Term Care

Follow your in-house disinfection procedures and only use the disinfectants and methods specified therein.

7.4 Service interval

WARNING!

Risk of injury or damage

Service must be performed only by a qualified technician.
 – Contact your Invacare provider for service.

Unless otherwise stated in local requirements, service must be performed at least every 12 months and according to the Service Inspection Checklist. When performing annual or regular service, all parts designed to carry load must be, as a minimum, tested with maximum load. All safety features must be checked according to ISO 10535:2006 Annex B.

7.4.1 Safety inspection checklist

A person who is suitably and properly qualified and well acquainted with the design, use and care of the lift should perform periodic inspections.

Date of Inspection:	Initials:
THE CASTOR BASE	
<input type="checkbox"/> Inspect for missing hardware. <input type="checkbox"/> Base opens/closes with ease. <input type="checkbox"/> Inspect castors and axle bolts for tightness. <input type="checkbox"/> Inspect castors for smooth swivel and roll. <input type="checkbox"/> Inspect and clear castors of debris. <input type="checkbox"/> Inspect pivot joints for wear.	
SLINGS AND HARDWARE	
<input type="checkbox"/> Check all sling attachments each time it is used to ensure proper connection and patient safety. <input type="checkbox"/> Inspect sling material for wear. <input type="checkbox"/> Inspect straps for wear. <input type="checkbox"/> Inspect stitching.	
THE HYDRAULIC PUMP/ELECTRIC ACTUATOR ASSEMBLY	
<input type="checkbox"/> Check for leakage. <input type="checkbox"/> Inspect hardware on mast, boom and base. <input type="checkbox"/> Check for wear or deterioration. If damaged, return to factory. <input type="checkbox"/> Cycle to ensure smooth quiet operation of the electric actuator. <input type="checkbox"/> Check for smooth operation of the pump handle for hydraulic pumps.	
THE BOOM	
<input type="checkbox"/> Check all hardware and hanger bar supports. <input type="checkbox"/> Inspect for bends or deflections. <input type="checkbox"/> Inspect bolted joints of boom for wear. <input type="checkbox"/> Inspect to ensure that the boom is centered between the base legs. <input type="checkbox"/> Check the mast pivot bolt. Ensure that the bolt is tightly secured. <input type="checkbox"/> Inspect pivot joints for wear.	

THE MAST
<input type="checkbox"/> Mast must be securely assembled to boom. <input type="checkbox"/> Inspect for bends or deflections. <input type="checkbox"/> Inspect pivot joints for wear.
THE SPREADER BAR
<input type="checkbox"/> Check the bolt/hooks for wear or damage. <input type="checkbox"/> Check sling hooks for wear or deflection. <input type="checkbox"/> Inspect pivot joints for wear.
LEG SPREADER HANDLE
<input type="checkbox"/> Ensure leg spreader handle operates smoothly.
CLEANING
<input type="checkbox"/> Whenever necessary.


7.5 Lubricating the Lift

WARNING!

Risk of Falling


Hydraulic oil or lubricant on the floor can cause a fall and injury.
 – Wipe excess lubricant from the lift after lubrication.
 – Use a facial tissue to wipe excess oil from around the hydraulic piston.
 – If excess oil is leaking from the hydraulic pump, contact a dealer or qualified technician for service.

Refer to the figure for lubrication points. Lubricate all pivot points with a light grease (waterproof auto lubricant). Wipe all excess lubricant from lift surface.

1. Spreader Bar
2. Boom Mounting Bracket
3. Boom/Mast Mount
4. Mast Mounting Bracket

7.6 Checking and tightening bolts

1. Check that the bolts are through the brackets and the locknuts are tight and secure.
2. If needed, do one or more of the following:
 - Tighten locknut and back-off the locknut 1/8 of a turn.
 - Replace the locknut.

7.7 Maintaining the Hydraulic Pump

(Atlante Hydraulic only)

WARNING!

Risk of Damage

The pump is sealed at the factory. To avoid voiding the warranty or causing product damage:

- If service is required, the pump unit **MUST** be returned to Invacare for repair.
- **DO NOT** attempt to open the pump or obtain local service. This will void the warranty and might result in damage and costly repair. Consult your dealer or Invacare for further information.

Risk of Falling

Hydraulic oil or lubricant on the floor can cause a fall and injury.

- Wipe excess lubricant from the lift after lubrication.
- Use a tissue to wipe excess oil from around the hydraulic piston.
- If excess oil is leaking from the hydraulic pump, contact a dealer or qualified technician for service.

All parts of the hydraulic pump are precision machined, then carefully assembled and tested to ensure reliable service. The pump assembly is completely enclosed and sealed with neoprene rings to prevent leakage of hydraulic oil. A small amount of oil (about a drop) will accumulate around the piston from time to time and should be removed with a tissue.

8 After Use

8.1 Transport and Storage

During transportation, or when the patient lift is not to be used for some time, the emergency stop button should be pushed in.

A disassembled lift can be located in the packaging box, for transportation and storage. See 4.4 *Disassembling the Lift*, page 10.

The patient lift must be stored at normal room temperature. If it is stored in a damp, cold or wet environment then the motor and other mounting parts may be prone to corrosion. See 10.4 *Environmental conditions*, page 25.

8.2 Disposal

WARNING!

Environmental Hazard

Device contains batteries.

This product may contain substances that could be harmful to the environment if disposed of in places (landfills) that are not appropriate according to legislation.

- DO NOT dispose of batteries in normal household waste.
- Batteries MUST be taken to a proper disposal site. The return is required by law and free of charge.
- Do only dispose discharged batteries.
- Cover terminals of lithium batteries prior to disposal.
- For information on the battery type see battery label or chapter 10 *Technical Data*, page 24.

Be environmentally responsible and recycle this product through your recycling facility at its end of life.

Disassemble the product and its components, so the different materials can be separated and recycled individually.

The disposal and recycling of used products and packaging must comply with the laws and regulations for waste handling in each country. Contact your local waste management company for information.

8.3 Reconditioning

This product is suitable for reuse. To recondition the product for a new user, carry out the following actions:

- Inspection
- Cleaning and disinfection

For detailed information, see 7 *Maintenance*, page 19.

Make sure that the user manual is handed over with the product.

If any damage or malfunction is detected, do not reuse the product.

9 Troubleshooting

9.1 Identifying and repairing faults

WARNING!

– Only personnel having received the necessary instruction or training by Invacare must perform service and maintenance on the product.


Symptoms	Faults	Solution
Patient lift feels loose.	Mast/Base joint loose.	Refer to Assembling the Mast to the Base in the Setup section.
	Tie - Rods are loose.	Contact your Dealer or Invacare representative.
Casters/Brakes noisy or stiff.	Fluff or debris in bearings.	Clean casters from fluff and debris.
Noisy or dry sound from pivots.	Needs lubrication.	Refer to Lubricating the Lift in the Maintenance section.
Electric actuator fails to lift.	Hand-control or actuator connector loose.	Connect hand control or actuator connector. Ensure connectors are seated properly and fully connected.
	Battery low.	Charge batteries. Refer to Charging the battery in the Usage section.
	RED emergency stop button pressed IN.	Rotate RED emergency stop button CLOCKWISE until it pops out.
	Battery not connected properly to control box.	Reconnect the battery to the control box. Refer to Charging the battery in the Usage section.
	The connecting terminals are damaged.	Replace the battery pack. Refer Charging the battery in the Usage section.
	Actuator needs service or load exceeds max. user weight.	Contact your Dealer or Invacare representative.
Unusual noise from actuator.	Actuator is worn or damaged or spindle is bent.	Contact your Dealer or Invacare representative.
Boom will not lower in uppermost position.	Boom requires a minimum weight load to lower from the uppermost position.	Pull down slightly on the boom.
Boom will not lower during a power retraction.	Shoulder bolt at the junction of the boom and mast may not be properly installed.	Refer to Checking and Tightening Bolts in the Maintenance section.
The control unit emits a beeping sound during lifting, and the motor stops.	Max. load is exceeded	Reduce the load (and the lift will function normally).
Manual/Hydraulic pump fails to lift when pumped.	Control valve not fully closed.	Close the control valve.
	Manual/Hydraulic pump in need of replacement	Contact your Dealer or Invacare representative.
Oil leaking from hydraulics.	Manual/Hydraulic pump in need of replacement.	Contact your Dealer or Invacare representative.

Contact your dealer or Invacare representative if the above does not solve your problems.

10 Technical Data

10.1 Dimensions and weight

* Forward direction

Dimensions [mm]	Kirilos	Atlante Electric	Atlante Hydraulic
Front / rear castor diameter	75	100	
Max. reach at 600 mm (a)	480		
Max. reach from base (b)	520		
Reach from base with legs spread to 700 mm (c)	260		
CSP* max. height / highest position (k)	1780	1820	
CSP* min. height / lowest position (l)	850	890	
Lifting range (m)	930		
Total width (open) centre to centre of castor (p)	960		
Internal width at maximum reach (q)	810		
Min. internal width (r)	580		
Min. free height (s)	30	60	
Height to upper edge of legs (t)	110	150	
Base length (u)	1170	1200	

Dimensions [mm]	Kirilos	Atlante Electric	Atlante Hydraulic
Max. internal base length (v)	870	890	
Total width (closed), external measure (w)	680	700	
Total width (closed), internal measure (x)	530	510	
Total width (open) external measure (y)	1030	1060	
Total width (open) internal measure (z)	900	870	
Turning radius	1230		
Min. height at Sling Hook-up	575	610	
Max. height at Sling Hook-up	1490	1530	

* CSP = Central Suspension Point

Weights [kg]	Kirilos	Atlante Electric	Atlante Hydraulic
Max. Safe Working Load (patient + spreader bar + sling)	140		
Total weight incl. spreader bar	33	37	
Weight, mast assembly, incl. battery, boom, spreader bar and actuator or pump	17	21	
Base section with legs	16		

10.2 Electrical system

	Kirilos	Atlante Electric
Voltage output	24 V DC, max. 250 VA	
Voltage supply	100 – 240 V AC, 50/60 Hz	
Maximum current input	400 mA	
Protection class (entire device)	IPX4*	
Insulation class	Class II equipment	
Type B applied part	Applied part complying with the specified requirements for protection against electrical shock according to IEC60601-1.	
Sound level	45 – 50 dB (A)	
Working ability	40 full lifts without battery charge with batteries at 50% of full capacity	
Intermittent (periodic motor operation)	10%, max. 2 min. / 18 min.	
Battery specifications	2 x 12 V / 2.9 Ah	
Battery type	Lead acid (service free, sealed)	

	Kirilos	Atlante Electric
Manual emergency lowering	Yes	
Electric emergency lowering / lifting	Yes	

*IPX4 - The system is protected against water splashed from any direction.

10.3 Hydraulic system

Technical specifications of hydraulic pump

	Atlante Hydraulic
Push rod max. stroke	300 mm
Push force (dynamic)	310 daN. at 95 Bar
Maximum static load	390 daN. at 120 Bar
Rod load / lever load ratio	20 / 1
Weight	9 Kg. (9.2 daN)

10.4 Environmental conditions

	Storage and transportation	Operation
Temperature	-10 °C to +50 °C	+5 °C to +40 °C

Relative humidity	20% to 90% at 30 °C, not condensing
Atmospheric pressure	700 hPa to 1060 hPa

Allow the product to reach operation temperature before usage:

- Warming up from minimum storage temperature might take more than 30 minutes.
- Cooling down from maximum storage temperature might take more than 5 minutes.

10.5 Materials

Component	Material
Base, legs, mast and boom	Steel, powder-coated
Spreader bar	Steel, powder-coated
Actuator housing, hand control, mast protector, castors and other plastic parts	Material according to marking (PA, PP, PE)
Boom connector, bolts and nuts	Steel, zinc-plated

All components of the product are either corrosion resistant or corrosion protected.

11 Electromagnetic compatibility (EMC)

11.1 General EMC information

Medical Electrical Equipment needs to be installed and used according to the EMC information in this manual.

This product has been tested and found to comply with EMC limits specified by IEC/EN 60601-1-2 for Class B equipment.

Portable and mobile RF communications equipment can affect the operation of this product.

Other devices may experience interference from even the low levels of electromagnetic emissions permitted by the above standard. To determine if the emission from this product is causing the interference, run and stop running this product. If the interference with the other device operation stops, then this product is causing the interference. In such rare cases, interference may be reduced or corrected by the following:

- Reposition, relocate, or increase the separation between the devices.

11.2 Electromagnetic emission

Guidance and manufacturer's declaration

This product is intended for use in the electromagnetic environment specified below. The customer or the user of this product should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group I	This product uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	This product is suitable for use in all establishments including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations /flicker emissions IEC 61000-3-3	Complies	

11.3 Electromagnetic Immunity

Guidance and manufacturer's declaration

This product is intended for use in the electromagnetic environment specified below. The customer or the user of this product should assure that it is used in such an environment.

Immunity test	Test / Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient / burst IEC 61000-4-4	± 2 kV for power supply lines; 100 kHz repetition frequency ± 1 kV for input / output lines; 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line to line ± 2 kV line to earth	Mains power quality should be that of a typical commercial or hospital environment.

Immunity test	Test / Compliance level	Electromagnetic environment – guidance
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 0% U_T for 0,5 cycle at 45° steps 0% U_T for 1 cycles 70% U_T for 25 / 30 cycles < 5% U_T for 250 / 300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of this product requires continued operation during power mains interruptions, it is recommended that the product is powered from an un-interruptible power supply or a battery. U_T is the a. c. mains voltage prior to application of the test level.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6	3 V 150 kHz to 80 Mhz 6 V in ISM & amateur radio bands	Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which this product is used exceeds the applicable RF compliance level above, this product should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating this product. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	10 V/m 80 Mhz to 2,7 GHz 385 MHz - 5785 MHz test specifications for immunity to RF wireless communication equipment refer to table 9 of IEC 60601-1-2:2014	Portable and mobile RF communications equipment should be used no closer than 30 cm to any part of this product including cables.

 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

11.4 Test specifications for immunity to RF wireless communications equipment

IEC 60601-1-2:2014 — Table 9

Test Frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380 - 390	TETRA 400	Pulse modulation ^{b)} 18 Hz	1.8	0.3	27
450	430 - 470	GMRS 460, FRS 460	FM ^{c)} ± 5 kHz deviation 1 kHz sine	2	0.3	28
710 745 788	704 - 787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	0.2	0.3	9
810 870 930	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	2	0.3	28

Test Frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum power (W)	Distance (m)	Immunity test level (V/m)
1720 1845 1970	1700 - 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	2	0.3	28
2450	2400 - 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	0.3	28
5240 5500 5785	5100 - 5800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0.2	0.3	9

 If necessary to achieve the immunity test level, the distance between the transmitting antenna and the Medical Electrical Equipment or System may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

^{a)} For some services, only the uplink frequencies are included.

^{b)} The carrier shall be modulated using a 50 % duty cycle square wave signal.

^{c)} As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

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Invacare Poirier SAS