# Invacare® LiNX

DLX-REM050, DLX-REM060, DLX-REM1xx, DLX-REM2xx

en	Joystick Gaiter Service Manual .....	1
de	Joystick-Manschette Servicehandbuch.....	4
es	Pieza de joystick Manual de servicio.....	6
fr	Soufflet d'étanchéité du joystick Manuel de maintenance.....	9
it	Cuffia per joystick Manuale per la manutenzione .....	12
pt	Fole do joystick Manual de assistência.....	14
sv	Joystickdamask Servicemanual .....	17

## 1 Before You Start

### Distribution

These instructions are intended for use by experienced service technicians only. Gaiters and associated parts are not supplied for retail sale or installation by end users.

### Warranty

This procedure is intended for the replacement of gaiters on units not covered by warranty. Undertaking this procedure on a product within warranty will void the warranty. All units under warranty should be returned to an authorised Dynamic Controls service centre for repair, refund or replacement under warranty.

### Liability

Dynamic Controls accepts no responsibility for any liability, whether expressed or implied, resulting from use of unauthorized parts or servicing that deviates from this procedure.

### Limitations

This procedure covers the replacement of gaiters only. If any secondary fault, including, but not limited to, the presence of internal contamination, is discovered during the process of the gaiter replacement, the unit must be withdrawn from service and either exchanged or returned to a Dynamic Controls service centre for further analysis and where appropriate repair.

### Disclaimer

The remote modules covered under this procedure are supplied with tamper evident warranty void seals fitted. If a remote module is received for gaiter replacement without the seal in an intact condition, Dynamic Controls recommends that the unit be returned 'as is' without

en

undertaking any service or repair. This recommendation is intended to safeguard the healthcare professional or service technician from the consequences of damage or defects resulting from user tampering, whether intentional or inadvertent. Following repairs, the technician undertaking gaiter replacement under this procedure may choose to apply their own tamper evident labels, to protect against subsequent opening or tampering of the unit post-service.

### ESD precautions

ESD precautions are not required for gaiter replacement in the REM050, REM060, REM1xx and REM2xx remote modules.

## 2 Equipment

### Parts List

 Different variants of the gaiter and gaiter retainer are available on the market, which differ from each other. Ensure that you use the enclosed part for repair only!

Quantity	Description
1	Gaiter (required)
1	DOWSIL (formerly Dow Corning) DC4 Electrical Insulating Compound (required)
1	Gaiter retainer (optional)

 DC4 electrical insulating compound can be sourced from: Element14 / Farnell — Order Code 537019

### Tools required

-  • 1 x TORX screwdriver bit, size T10  
 • 1 x torque driver capable of taking the T10 TORX bit and reaching at least 80 cNm torque  
 • Cotton buds / swabs

DEALER: Keep this manual.  
The procedures in this manual MUST be performed by a qualified technician.

Yes, you can.®

### 3 Initial Checks

Before removing, the gaiter should be checked carefully to confirm the extent of wear or damage. If the gaiter is damaged to an extent that may have allowed foreign objects or significant ingress to pass into the joystick housing, the remote module should be returned for thorough internal assessment by a Dynamic Controls service centre.

If the remote module shows any signs of significant mechanical damage or impact, the unit should be returned to a Dynamic Controls service centre for assessment. Significant damage would typically include loss of function or openings in the unit chassis or speed dial.

### 4 Disassembly

**WARNING!**

– Do not disassemble beyond the point described in this procedure. Further disassembly may result in product calibration moving out of specification, compromising the product’s safe operation.

1. Remove the joystick knob from the remote module by pulling it directly away from the remote module's body. The knob may be held tightly on the shaft, so some force may be required — **DO NOT** twist the knob while it is on the shaft.

**WARNING!**

– Twisting the knob when it is still attached to the shaft may result in damage to the joystick mechanism.

2. Using the T10 driver, remove the four screws, highlighted in images below, from the underside of the remote module — the screws can be removed in any order. For REM060, REM1xx and REM2xx modules, one screw is located underneath the warranty label. **DO NOT** remove any other screws.

Fig. 4-1 REM050 gaiter retainer screw locations

Fig. 4-2 REM060, REM1xx and REM2xx gaiter retainer screw locations. Fourth screw **A** under label.

3. Pull the gaiter and gaiter retainer away from the top case and joystick shaft and then separate the gaiter from gaiter retainer.


Fig. 4-3 REM050 gaiter and gaiter retainer


Fig. 4-4 REM060, REM1xx and REM2xx gaiter and gaiter retainer

### 5 Component Inspection

**Contamination**

With the gaiter and retainer removed, inspect the visible joystick mechanism for signs of ingress, including signs of liquid residue and solid debris.

Deflect the joystick shaft in all directions so that the inside of the joystick mechanism can also be checked for signs of debris or ingress. If there is any sign of contamination, the complete remote module should be replaced or, where appropriate, returned for service by a Dynamic Controls service centre.

**WARNING!**

– Unauthorised cleaning of the joystick mechanism or application of unsuitable lubricants may adversely affect the performance of the joystick mechanism, creating a safety risk for the user.


**Spring Retainer And Shaft**

Fig. 5-1 Spring retainer **A**, spring retainer pin **B** and shaft

While the gaiter is removed, check the condition of the spring retainer and shaft. Check for signs of physical damage, including deformed shafts and cracks in the spring retainer. Ensure the spring retainer's locating pin is secure and not protruding from the retainer, as shown in 5-1 *Spring retainer **A**, spring retainer pin **B** and shaft, page 3*. Note that the mould line on the spring retainer is normal and not a sign of damage.

If any part is damaged, loose or missing, the complete remote module should be replaced or, where appropriate, returned for service by a Dynamic Controls service centre.

**Gaiter Retainer**

Inspect the gaiter retainer for any signs of damage, including cracking of the posts. If any damage is visible, the retainer should be replaced.

**6 Reassembly**

1.

Fig. 6-1 Gaiter Pushed Through The Gaiter Retainer

Push the gaiter through the gaiter retainer.

2.


Fig. 6-2 Middle Layer


Fig. 6-3 Lower Layer


Fig. 6-4 Cross Section (Middle And Lower Layers)

Using a cotton bud, apply a thin coating<sup>1</sup> of DC4 Electrical Insulating Compound to the middle and lower layers of the gaiter, as highlighted with the grid in figures above. Ensure that the grease is smeared evenly all the way around the highlighted areas and not placed in just one spot.

3.


Fig. 6-5 Push the gaiter and gaiter retainer onto the top case assembly


Fig. 6-6 Check gaiter sits correctly over spring retainer **A**

Push the gaiter and gaiter retainer onto the top case assembly, whilst rotating the gaiter within the retainer to ensure that it can move freely. Ensure that the gaiter sits down fully and does not have any folds or wrinkles and that it sits correctly on the spring retainer as shown in 6-6 *Check gaiter sits correctly over spring retainer **A**, page 3*.

1. The amount of insulating compound that you use is not critical. Use just enough so that the surfaces, as indicated in the images, are lightly covered.

4. Replace the four screws used to secure the gaiter retainer into the bottom case and tighten to 80 cNm — the screws can be replaced in any order. For security, dealers may wish to attach their own anti-tamper labels so that they can determine if the unit is subsequently opened.
- 5.


Fig. 6-7

Replace the joystick knob. If necessary, rotate the knob to align with the notches **A** and then push the knob firmly towards the body of the remote module— when fitted correctly, the knob should not be able to be rotated about the shaft.

## 7 Testing

Visually inspect the gaiter's position and ensure the gaiter retainer sits flat with the top case: the gaiter must sit down fully, with no folds or wrinkles in it.

Test the joystick movement by moving the joystick in a circular motion about its full extent. It must move freely in all directions and consistently return to the centre, neutral position when released from any position.

Connect to a system. Check for any out-of-neutral at power up or other faults.

Use the remote module to drive a wheelchair, ensuring that the wheelchair drives and stops correctly, the joystick returns to neutral and that, after a short while, the park brakes engage (usually indicated by an audible click).

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Invacare Ireland Ltd,
Invacare AB
Invacare Deutschland GmbH