**C300T40
C300T44
C333T40
C333T43
C333T45
(+ C333-00001 et
C333-00002)
L880
L883
(+ L881)**

F	1.	Guide utilisateur
GB	2.	User's manual
D	3.	Gebrauchsanweisung
E	4.	Manual del usuario
I	5.	Manuale dell' utente
NL	6.	Gebruiksaanwijzing
P	7.	Guia do utilizador
S	8.	Bruksanvisning
DK	9.	Brugsanvisning
FIN	10.	Käyttäjän käsikirja
N	11.	Bruksanvisning

I. FRANCAIS

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
2. ENGLISH

You have just purchased an INVACARE® product. Our products are designed and manufactured to the most rigorous quality and safety standards and conform to the European directives (MDD 93/42/EEC). They are guaranteed for one year. Certain products may have a longer guarantee; please check with your dealer.

The following operating and maintenance instructions must be observed:

- Check that the product has been delivered in perfect condition before you use it.
- If you have any doubts about the restrictions on use, the adjustment or maintenance of the product, please contact your dealer.
- Do not exceed the maximum operating load  or the maximum load  shown on the product label.
- Always ensure that the product is suitable for your requirements; in case of doubt, check with a health professional.
- Worn, damaged or missing components must be replaced immediately by an approved repairer.
- Do not leave or store the product out of doors.
- Periodically check the stability and rigidity of the product and, where applicable, that the screws are correctly tightened and the rubber end caps in good condition.
- For reasons of hygiene, the same person should always use the mattress and pillows.
- As a primary precaution against bedsores, change the body's support points (as often as possible for the pillows and every three hours for the mattress) to avoid compressing blood vessels for too long at the bony protuberances and to assist circulation.

C300T40 - C300T44 - C333T40 - C333T43 - C333T45 PILLOWS AND C333-00001 - C333-00002 PILLOWCASES

- The pillows are designed for users weighing between 40 and 100 kg.
- For health reasons, pillows should be replaced after one year's use.
- C333: Slip the pillow into the pillowcase with the non-slip side uppermost.

Maintenance

Wash the pillowcase in pH neutral liquid detergent and disinfect with diluted, denatured alcohol (water/alcohol ratio 3:1). Tumble dry at 50°C max. Store flat between 0 and 40°C.

L880 AND L883 CELLULAR MATTRESSES L880 AND L881 UNDERSHEETS

- These mattresses are designed for users weighing between 40 and 100 kg.
- To maintain their efficiency the overmattresses should be covered with a flat undersheet and top sheet.

L883 multidensity mattress

This mattress is composed of a green base and 9 removable modules positioned in the following way:

- Beige: head (soft and comfortable main contact points).
- Yellow: very firm foam. The main contact points are positioned at the body's support points, i.e. the buttocks and heels (to protect against bedsores).
- Pink: legs and back (medium hard).

Maintenance

Clean the mattress with a clean cloth and soapy water, then dry with an absorbent terry towelling cloth. The hypoallergenic cover can be cleaned with normal washing powder and disinfected by boiling at temperatures up to 90° or using chlorinated or formalin-based products; do not use solvents. The cover can also be cleaned and disinfected by soaking for 15 minutes in a BACTINYL or similar solution. Store the mattress in a clean dry place.

USER'S MAXIMUM WEIGHT

Please see the table at the end of this document.

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